



For research use only

ISO9001

His•Bind Magnetic Agarose Resin (Ni-IDA)

Product Name	Qty	Cat. No.	Remarks
His•Bind Magnetic Agarose Resin (Ni-IDA)	2 ml	EBE-1032	pre-charged with Ni ion, 75% slurry

Description

Elpis Biotech's His•Bind agarose resin is designed for rapid purification of recombinant proteins containing His•Tag sequence by IMAC. The His•Tag sequence binds to Ni²⁺ cations, which are immobilized on the His•Bind agarose resin. Para-magnetic particles are embedded onto agarose bead enables to purify in a batch style high-throughput purification using a magnet.

After unbound proteins are washed away, the target protein is recovered by elution with imidazole containing wash buffer. The His•Bind agarose resin can be regenerated and reused several times.

Features

- **Matrix** Agarose 4B, with magnetite, cross-linked
- **Activation** Epoxy
- **Ligand** Iminodiacetic acid (Ni ion pre-charged)
- **Linker Size** 12-atom space linker
- **Binding Capacity** 5-10 mg polyhistidine tagged proteins/ml resin
- **Bead Size** 10-25 μm
- **pH Stability** pH 2 to 14
- **Storage** 75% suspension in PBS with 20% ethanol

Storage

Store at 4°C, do not freeze

Reagent Compatibility

β -Mercaptoethanol	< 10 mM
CHAPS	< 1%
Ethanol	30%
Ethylene glycol	30%
HEPES	50 mM
Glycerol	20%
Guanidine HCl	6 M
Imidazole	500 mM at pH 7.0–8.0, for elution
KCl	500 mM
MES	20 mM
MOPS	50 mM
NaCl	1.0 M
NP-40	1%
SDS	< 1%
TRIS	50 mM
Triton-X 100	<1%
Urea	8 M

• Reagent compatibility shows the maximum concentrations of each reagent tested. Higher levels may be acceptable, but they should be tested before use. Note that some of these reagents may partially or completely denature your protein. EDTA or EGTA can not be used completely.

Protocol

Resin preparation

1. Prepare 1x Binding buffer, 1x Washing buffer, and 1x Elution buffer before purification (see below for buffer composition).
2. Gently mix the tube of His•Bind magnetic agarose resin by inversion until completely suspended. Using a wide-mouth pipette, transfer the desired amount of agarose resin slurry to a column.
3. Wash the resin with 3 volumes 1x binding buffer.

** Caution : The resins are supplied as pre-charged form with Ni ion.

Small Scale Purification

1. Resuspend the cells in 1X Bind buffer at 2–5 ml per gram wet weight. The amount of cells required depends on the expression level of the His•Tag fusion protein and the expression system used. The binding capacity of His•Bind magnetic agarose resins is protein-dependent and normally lies between 5–10 mg/ml. Lysis/binding buffer contains 5 mM imidazole to minimize binding of untagged, contaminating proteins and increase purity with fewer wash steps. If the His•Tag fusion protein does not bind under these conditions, the amount of imidazole should be reduced to 1–5 mM. With His•Tag fusion proteins exhibiting high binding affinities, the imidazole concentration can be increased to 20 mM.
2. Add lysozyme to 1 mg/ml and incubate on ice for 30 min
3. Sonicate on ice. Use six 10 sec bursts at 200–300 W with a 10 sec cooling period between each burst. Use a sonicator equipped with a microtip
4. (Optional) If the lysate is very viscous, add RNase A (10 µg/ml) and DNase I (5 µg/ml) and incubate on ice for 10–15 min.
5. Centrifuge lysate at 13,000rpm for 20–30 min at 4°C to pellet the cellular debris. Save supernatant and store –20°C.
6. Add 100-500 µl of cleared lysate to the 20 µl His•Bind magnetic agarose resin and mix gently by shaking (200rpm on a rotary shaker) at 4°C for 60 min.
7. Centrifuge (optionally by magnetic separator) lysate- His•Bind magnetic agarose mixture at 3,000rpm for 30 sec at 4 °C.
8. Remove supernatants. Save flow-through for SDS-PAGE analysis.
9. Wash twice with 500 µl 1x Wash buffer by repeat of mixing and centrifugation at 3,000rpm for 30 sec at 4 °C, collect wash fractions. Save wash fractions for SDS-PAGE analysis.
10. Elute the protein with 4 × 20 µl 1x Elution buffer. Collect the eluate in four tubes and analyze fractions by SDS-PAGE. Alternatively, 10 -100 µl of Strip buffer may also be used to remove the protein by stripping the Ni²⁺ from the column.

Buffers Used

1. Buffers for purification under native conditions

1X Bind buffer (use for lysis and binding): 0.5 M NaCl, 5 mM imidazole, 20 mM Tris-HCl, pH 7.9

1X Wash buffer: 0.5 M NaCl, 60 mM imidazole, 20 mM Tris-HCl, pH 7.9

1X Elute buffer: 0.25 M NaCl, 0.5 M imidazole, 10 mM Tris-HCl, pH 7.9

1x Stripping buffer: 0.5 M NaCl, 100 mM EDTA, 20 mM Tris-HCl, pH 7.9

1x Charging buffer: 50 mM NiSO₄

2. Buffers for purification under denaturing conditions

Denaturing Lysis/Bind buffer: 6 M Guanidine-HCl, 0.1 M NaH₂PO₄, 0.01 M Tris-Cl, pH 8.0
(optionally 8 M urea, 0.1 M NaH₂PO₄, 0.01M Tris-Cl, pH 8.0)

Denaturing Wash buffer: 8 M urea, 0.1 M NaH₂PO₄, 0.01 M Tris-Cl, pH 6.3

Denaturing Elution buffers: 8 M urea; 0.1 M NaH₂PO₄; 0.01 M Tris-Cl pH 5.9
(optionally 8 M urea, 0.1 M NaH₂PO₄, 0.01 M Tris-Cl, pH 4.5)