



For research use only

ISO9001

## MBP•Bind Magnetic Agarose Resin

Product Name	Qty	Cat. No.	Remarks
MBP•Bind Magnetic Agarose Resin	2 ml	EBE-1072	75% slurry

### Description

Elpis Biotech's MBP•Bind magnetic agarose resin is designed for ultra rapid purification of recombinant proteins containing MBP (maltose binding protein) tag sequence. The MBP•Tag sequence binds to amylose, which is immobilized on the MBP•Bind magnetic agarose resin. The high degree of substitution of amylose to magnetic agarose resin ensures a high binding capacity with yields of MBP fusion proteins of 5-10 mg/ml settled resin. Para-magnetic particles are embedded onto agarose bead enables to purify in a batch style high-throughput purification using a magnet. The MBP•Bind magnetic agarose resin is suitable for small scale purification rather than large purification.

After unbound proteins are washed away, the target protein is recovered by elution with maltose containing elution buffer. The MBP•Bind resin can be regenerated and reused many times.

### Features

- **Matrix** Agarose 4B with magnetite, cross-linked
- **Activation** Epoxy
- **Ligand** Hydroxyl, amylose
- **Linker Size** 12-atom space linker
- **Bead Size** 10-25  $\mu\text{m}$
- **Binding Capacity** 5-10 mg MBP•Tagged proteins/ml resin
- **pH Stability** pH 2 to 14
- **Storage** 75% slurry in PBS with 20% ethanol

### Storage Condition

Store at 4°C, do not freeze

### Resin preparation

1. Prepare 1x Bind/washing buffer, and 1x Elution buffer before purification (see below for buffer composition).
2. Gently mix the tube of MBP•Bind magnetic agarose resin by inversion until completely suspended. Using a wide-mouth pipette, transfer the desired amount of magnetic agarose resin slurry to a tube.
3. Wash the bead with 3 volumes of 1x Bind/washing buffer by centrifugal or magnetic force.

### Small Scale Purification

1. Add pre-washed magnetic beads to the prepared extract (50  $\mu\text{l}$  to 1 ml) in 1x Bind/washing buffer. Allow binding MBP•tagged proteins to MBP•Bind magnetic agarose resin for 5-10 min at room temperature or at 4°C with gentle agitation (10 –20  $\mu\text{l}$  bead volume is sufficient to purify 50-100  $\mu\text{g}$  MBP tagged proteins, but the optimized amount of bead volume must be determined empirically).
2. Wash the bead three times with 200  $\mu\text{l}$  1X Bind/washing buffer by centrifugal or magnetic force.
3. Elute the bound protein with 10 – 100  $\mu\text{l}$  of 1X Elute buffer.

### Buffers Used

Bind/washing buffer : 20 mM Tris-HCl, pH 7.4 , 200 mM NaCl, 1 mM EDTA, 10 mM  $\beta$ -mercaptoethanol

Elution buffer : Bind/washing buffer + 10 mM maltose