



For research use only

ISO9001

100bp DNA Ladder Marker, Ready-to-use

Product	Conc.	Cat. No.	Remarks
100bp DNA Ladder Marker, Ready-to-use	1 ml (33 µg/ml)	EBM-1001	Ready-to-use

Description

The 100bp DNA ladder marker is a mixture of specially designed double-stranded DNA fragments for determining the exact size of PCR products and engineered DNA fragments. The 100bp DNA ladder marker consists of 11 DNA fragments ranging in size from 100 to 1,000 bp in 100 bp increment, and additional 1,500 bp fragment. For easy size reference on the gel electrophoresis, the 500 bp and 1,000 bp are two to three times more brighter than the other bands. The 100bp DNA ladder marker is supplied in a ready-to-use format. This ladder marker can be stained with ethidium bromide or any other known DNA staining methods.

Storage Buffer

- Marker DNA : 33 µg in 1 ml of 10 mM Tris-HCl, pH8.0, 1 mM EDTA, 5% Glycerol, 0.005% Bromophenol Blue, and 0.005% Xylene Cyanol

Recommended Storage Condition

- -20°C for 2 year
- 4°C for 6 months
- Room temperature (20-25°C) for 2 months


Usage Information

- Concentration : 330 ng/10 µl (33 µg/ml)
- Recommended loading : 5-10 µl (100-200 lanes, ready-to-use)
- Range : 100 – 1,500 bp
- Number of bands : 11

Cautions

- Always use the fresh tip to take out marker solution.
(If you do not, trace amount of contaminated DNases from buffer tank may degrade marker DNA rapidly)
- Don't boil the product.
- Use appropriate % of gels for separation of 100 to 1,500 bp sizes
(1 to 3% agarose gel is recommended)
- Confirm that the concentration of DNA staining dye is optimal before use.
(Breakage or suboptimal concentration of ethidium bromide in gel is a main cause of low estimation of marker concentration or your DNA. 5 ng of DNA should be seen in normal condition)
- Loading volume and concentration should be optimized by gel size, well size, and running length.
- Low sized DNA bands can be gradually disappeared as running is progressing.
(This is because some DNA is getting out from gel to buffer during horizontal electrophoresis, not because the DNA concentration is incorrect. This will be the same for your DNA)


100bp DNA Ladder Marker



10 µl/330 ng/lane ;

2% agarose in 0.5x TBE, stained with ethidium bromide

Migration Patterns in Different % of Agarose Gels



Recommended Gel Percentages for Separation of Linear DNA

Agarose Gel, %	Range of Separation, bp	Polyacrylamide Gel, %	Range of Separation, bp
0.5	1,000 - 30,000	3.5	100 - 1,000
0.7	800 - 12,000	5	80 - 500
1	500 - 10,000	8	60 - 400
1.2	400 - 7,000	12	40 - 200
1.4	200 - 4,000	20	5 - 100
2	50 - 2,000		

DNA Size Migration with Sample Loading Dyes

Agarose Concentration, %	Xylene cyanol FF	Bromophenol blue	Orange G
0.7 - 1.7	~4000 bp	~300 bp	~50 bp
2.5 - 3.0	~800 bp	~100 bp	~30 bp

Composition of Gel Electrophoresis Buffers

Buffer	Working Concentration	Stock Concentration (per Liter)		
Tris-acetate (TAE)	1x	20 mM Tris-acetate	Tris base	96.9 g
		1 mM EDTA	Glacial acetic acid	22.84 ml
			0.5 M EDTA (pH8.0)	40 ml
Tris-borate (TBE)	0.5x	45 mM Tris-borate	Tris base	108 g
		1 mM EDTA	Boric acid	55 g
			0.5 M EDTA (pH8.0)	40 ml