



For research use only

ISO9001

HiPi Plus DNA Polymerase

Product Description	Quantity	Cat. No.	Remarks
HiPi Plus DNA Polymerase	500 unit	EBT-1007	5 unit/ μ l
	500 unit	EBT-1010	with 1 ml dNTP mix

Description

HiPi Plus DNA Polymerase is suitable for a high fidelity amplification of DNA fragments. This enzyme is designed for a reliable amplification of long (up to 10 kbp for lambda DNA) and complex targets with a robust yield and high specificity. HiPi Plus DNA Polymerase generates a mixture of PCR products with blunt end and 3'-dA overhangs.

Storage Buffer

5 unit/ μ l in 20 mM Tris-HCl, pH8.0, 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Tween 20, and 50% Glycerol.

Application(s)

High fidelity PCR, RT-PCR, genomic PCR, long PCR.

Unit Definition

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nmole of dNTP into an acid-insoluble form in 30 min at 72°C. The reaction conditions are : 25 mM TAPS, pH9.3, 50 mM KCl, 2 mM MgCl₂, 1 mM β -mercaptoethanol, 200 μ M each dNTPs, 100 μ M [α -³²P]dCTP, 12.5 μ g activated calf thymus DNA in a total volume of 50 μ l.

10x Reaction Buffer

500 mM Tris-HCl, pH9.0, 160 mM (NH₄)₂SO₄, 35 mM MgCl₂, 1% Tween 20, 1 mg/ml BSA.

QC tests

Activity, SDS-PAGE purity, performance tests.

Storage Condition

Store at -20°C.

Standard Protocol

1. Prepare 50 μ l PCR solution as follows:

PCR grade distilled water :	- μ l
10x reaction buffer :	5 μ l
10 mM dNTP mix (2.5 mM each) :	4 μ l
Primer (10 pmol/ μ l) :	1 μ l each
Template :	1-10 ng
HiPi Plus DNA Polymerase :	0.5 μ l (2.5 unit)

Adjust final vol. to 50 μ l with PCR grade distilled water

* Always, HiPi Plus DNA Polymerase should be added last to the mixture

* If you are amplifying a fragment larger than 2 kbp and with high GC content, add 5 μ l of 10x Q buffer into the PCR mix. It will greatly improve reaction specificity. Q buffer is helpful for GC rich template and long template.

2. Set PCR cycling as follows :

Initial denature at 95°C : 3 min

		< 1 kbp	1-10 kbp
Denature	95°C	30 sec	30 sec
Anneal	T _m -4°C	30 sec	30 sec
Extend	72°C	45 sec	30-60 sec/kbp

* 25-45 PCR cycles

* You can also use two step cycle for > 5 kbp amplification (denaturation at 95°C and annealing/extension at 68°C)

Trouble-Shooting

1. No products

- Confirm your template is intact : Try another reaction with a result assured primer pair and templates
- Be sure all the component are correctly added and working well : Sometimes low graded dNTP may inhibit the reaction, and degraded primers can result in low sized PCR fragments

2. Smear bands or smeared background

- Reduce template concentration : High concentration of template can lead to smearing of PCR products. Generally, 1-10 ng of plasmid DNA is working well
- Increase annealing temperature

3. Non-specific bands

- Increase annealing temperature
- Consider using PCR additives, like 1-2% DMSO and 0.5-1x Q buffer
- Confirm specificity of your primers

4. Low yield

- Increase enzyme concentration in the reaction
- Increase PCR cycle number
- Be sure appropriate concentration of your template is added