



For research use only

ISO9001

## Pfu DNA Polymerase

| Product Description | Quantity | Cat. No. | Remarks            |
|---------------------|----------|----------|--------------------|
| Pfu DNA Polymerase  | 500 unit | EBT-1011 | 5 unit/ $\mu$ l    |
|                     | 500 unit | EBT-1012 | with 1 ml dNTP mix |

### Description

Recombinant Pfu DNA polymerase is purified from an *E.coli* strain carrying a plasmid with the cloned gene encoding *Pyrococcus furiosus* DNA polymerase. The enzyme catalyzes the incorporation of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg<sup>2+</sup> at 70–80°C. Pfu DNA Polymerase exhibits 3'→5' exonuclease (proof-reading) activity, but has no detectable 5'→3' exonuclease activity.

The amplified products by Pfu DNA polymerase can be used for a gene cloning with decreased error rate, and for a site-specific mutagenesis. Pfu DNA Polymerase, like any other polymerases showing proof-reading activity, generates PCR products with blunt end.

Pfu DNA Polymerase is recommended for an amplification of DNA fragment smaller than 7 kbp. Pfu DNA Polymerase is provided with 10x optimized reaction buffer.

### Storage Buffer

5 unit/ $\mu$ l in 50 mM Tris-HCl, pH8.2, 0.1 mM EDTA, 1 mM DTT, 50% Glycerol.

### Unit Definition

One unit of enzyme catalyzes the incorporation of 10 nanomoles of deoxynucleotides into a polynucleotide fraction in 30 min at 72°C.

### 10x Reaction Buffer

200 mM Tris-HCl, pH9.0, 100 mM KCl, 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20 mM MgSO<sub>4</sub>, 1% Triton X-100, 1 mg/ml BSA.

### QC tests

Activity, SDS-PAGE purity, performance tests, genomic DNA contamination test, confirmation test for the absence of endo and exonucleases.

### Storage Condition

Store at -20°C.

### Standard Protocol

1. Prepare 50  $\mu$ l PCR solution as follows:

|                                |                        |
|--------------------------------|------------------------|
| PCR grade distilled water :    | - $\mu$ l              |
| 10x reaction buffer :          | 5 $\mu$ l              |
| 10 mM dNTP mix (2.5 mM each) : | 4 $\mu$ l              |
| Primer (10 pmol/ $\mu$ l) :    | 1 $\mu$ l each         |
| Template :                     | 1–10 ng                |
| Pfu DNA Polymerase :           | 0.5 $\mu$ l (2.5 unit) |

Adjust final vol. to 50  $\mu$ l with PCR grade distilled water

\*Note : Always, Pfu DNA polymerase should be added last to the mixture

2. Set PCR cycling as follows :

Initial denature at 95°C : 3 min

|          |                     | <1 kbp | 1–7 kpb         |
|----------|---------------------|--------|-----------------|
| Denature | 95°C                | 30 sec | 30 sec          |
| Anneal   | T <sub>m</sub> -4°C | 30 sec | 30 sec          |
| Extend   | 72°C                | 45 sec | 30–60 sec/1 kbp |

25–40 PCR cycles

### Trouble-Shooting

- No products
  - Confirm your template is intact : Try another reaction with a result assured primer pair and templates
  - Be sure all the component are correctly added and working well : Sometimes low graded dNTP may inhibit the reaction, and degraded primers can result in low sized PCR fragments.
- Smear bands or smeared background
  - Reduce template concentration : High concentration of template can lead to smearing of PCR products. Generally, 1–10 ng of plasmid DNA is working well
  - Increase annealing temperature
- Non-specific bands
  - Increase annealing temperature
  - Consider using PCR additives, like 1–2% DMSO and 0.5–1x Q buffer
  - Confirm specificity of your primers
- Low yield
  - Increase enzyme concentration in the reaction
  - Increase PCR cycle number
  - Be sure appropriate concentration of your template is added
- Mutation is found
  - Increase initial template concentration
  - Reduce PCR cycle number
  - Reduce dNTP concentration added in PCR mix