



For research use only

ISO9001

rTaq DNA Polymerase

Product Description	Quantity	Cat. No.	Remarks
rTaq DNA Polymerase	500 unit	EBT-1013	5 unit/ μ l
	500 unit	EBT-1014	with 1 ml dNTP mix

Description

Recombinant rTaq DNA polymerase is overexpressed and purified from an *E.coli* strain containing the gene of *Thermus aquaticus* DNA polymerase. rTaq DNA polymerase is modified to increase thermostability and performance in conventional PCR process.

Storage Buffer

5 unit/ μ l in 20 mM Tris-HCl, pH8.0, 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Tween 20, and 50% Glycerol.

Unit Definition

One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nmole of dNTP into an acid-insoluble form in 30 min at 72°C. The reaction conditions are : 25 mM TAPS, pH9.3, 50 mM KCl, 2 mM MgCl₂, 1 mM β -mercaptoethanol, 200 μ M each dNTPs, 100 μ M [α -³²P]dCTP, 12.5 μ g activated calf thymus DNA in a total volume of 50 μ l.

10x PCR Reaction Buffer w/ MgCl₂

500 mM Tris-HCl, pH9.0, 160 mM (NH₄)₂SO₄, 25 mM MgCl₂, 1% Tween 20, 1 mg/ml BSA.

QC tests

Activity, SDS-PAGE purity, performance tests.

Storage Condition

Store at -20°C.

Standard Protocol

1. Prepare 50 μ l PCR solution as follows:

PCR grade distilled water :	- μ l
10x reaction buffer :	5 μ l
10 mM dNTP mix (2.5 mM each) :	4 μ l
Primer (10 pmol/ μ l) :	1 μ l each
Template :	1-10 ng
rTaq DNA Polymerase :	0.5 μ l (2.5 unit)

Adjust final vol. to 50 μ l with PCR grade distilled water

*Note : Always, rTaq DNA polymerase should be added last to the mixture

2. Set PCR cycling as follows :

Initial denature at 95°C : 3 min

		<1 kbp	1-4 kbp
Denature	95°C	30 sec	30 sec
Anneal	T _m -4°C	30 sec	30 sec
Extend	72°C	45 sec	30-60 sec/1 kbp

25-40 PCR cycles

Trouble-Shooting

1. No products

- Confirm your template is intact : Try another reaction with a result assured primer pair and templates
- Be sure all the component are correctly added and working well : Sometimes low graded dNTP may inhibit the reaction, and degraded primers can result in low sized PCR fragments

2. Smear bands or smeared background

- Reduce template concentration : High concentration of template can lead to smearing of PCR products. Generally, 1-10 ng of plasmid DNA is working well
- Increase annealing temperature

3. Non-specific bands

- Increase annealing temperature
- Consider using PCR additives, like 1-2% DMSO and 0.5-1x Q buffer
- Confirm specificity of your primers

4. Low yield

- Increase enzyme concentration in the reaction
- Increase PCR cycle number
- Be sure appropriate concentration of your template is added