




DNA Polymerase I (*E.coli*)

Product	Quantity	Cat. No.	Remarks
DNA Polymerase I (<i>E.coli</i>)	200 unit	EBT-1022	5 unit/ μ l

Description

DNA Polymerase I from *E.coli* catalyzes the template-directed polymerization of nucleotides into duplex DNA in a 5'→3' direction. DNA Polymerase I possesses a 3'→5' exonuclease activity or "proof-reading" function, which lowers the error rate during DNA replication, and also contains a 5'→3' exonuclease activity, which enables the enzyme to replace nucleotides in the growing strand of DNA by nick translation. The enzyme, purified from recombinant *E.coli*, is capable of catalyzing *de novo* synthesis of synthetic homopolymers and provides a convenient method for the preparation of a variety of defined DNA substrates.



Concentration & Storage Condition

5 unit/ μ l. Store at -20°C.

Storage Buffer

50 mM Tris-HCl, pH 7.5, 1 mM DTT, 0.1 mM EDTA, 50% (v/v) glycerol.

10x Reaction Buffer

500 mM Tris-HCl, pH 7.2, 100 mM MgSO₄, 1 mM DTT.

Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmole of dNTP into TCA-insoluble material in 30 min at 37°C. The reaction conditions are: 67mM potassium phosphate, pH 7.5, 6.7 mM MgCl₂, 1 mM DTT, 50 μ g/ml activated calf thymus DNA and 33 mM dATP, dCTP, dGTP and dTTP (a mix of unlabeled and [³H]dTTP).

QC Tests

Activity, exo and endonuclease activity test, SDS-PAGE purity, performance tests.

Usage Information

1. Fill-In of 5'-Overhang to form Blunt ends

- Add the following components to the microcentrifuge tube :

DNA (1–5 μ g digested DNA containing 5'-overhangs)	x μ l
10x Reaction Buffer	2 μ l
1mM dNTP mix (0.25mM each)	5 μ l
DNA Polymerase I (2.5 units)	0.5 μ l
Nuclease-Free Water to final volume	20 μ l
- Incubate at 37°C for 1 hour.
- Heat at 75°C for 10 minutes to inactivate the enzyme.

2. Removal of 3'-Overhang to form Blunt ends

The 3'-overhang is first removed by the exonuclease activity of DNA Polymerase I.

- Add the following components to the microcentrifuge tube :

DNA (1–5 μ g digested DNA containing 5'-overhangs)	x μ l
10x Reaction Buffer	2 μ l
DNA Polymerase I (2.5 units)	0.5 μ l
Nuclease-Free Water to final volume	20 μ l
- Incubate at 37°C for 10 min.
- Add 5 μ l of the 1mM dNTP mixture to the DNA.
- Incubate at 12°C for 1 hour.
- Heat at 75°C for 10 min to inactivate the enzyme.

3. Nick Translation

This reaction may be scaled between 10-100 μ l volume. Nucleotide mix prepared by mixing equal volumes of the 3 unlabeled 300 μ M nucleotides chosen minus the nucleotide selected as label.

- Set up the following reaction in a microcentrifuge tube :

Nucleotide mix	10 μ l
Nick translation 10X buffer	5 μ l
Sample DNA (at 0.2 μ g/ μ l)	5 μ l
[α - ³² P]dCTP (400Ci/mmol at 10mCi/ml)	7 μ l
DNA Polymerase I/DNase I mix.	5 μ l
Nuclease-Free Water to final volume	50 μ l
- Incubate at 37°C for 1 hour.
- Add 5 μ l stop solution (0.25M EDTA (pH 8.0)).

Nick translation 10X buffer

500mM Tris-HCl (pH 7.2)
100mM MgSO₄
1.0mM DTT

DNA Polymerase I/DNase I mix

50% glycerol
50mM Tris-HCl (pH 7.2)
10mM MgSO₄
0.1mM DTT
0.5mg/ml nuclease-free BSA
1,000u/ml DNA Polymerase I
3u/ml RNase-Free DNase (DNase I)

Prepare the buffer solution and then add the DNA Polymerase I and RNase-Free DNase.