



For research use only

ISO9001

10x Tris-Glycine Running Buffer

Product	Volume	Cat. No.
10x Tris-Glycine Running Buffer	1 L (10x)	EBA-1042

Description.

Elpis Biotech's Tris-Glycine Running Buffer (10x) is suitable for use as an electrophoresis buffer for Tris-Glycine gel to separate proteins and for western blotting. Premixed electrophoresis buffer eliminate variables which can distort data and give unreliable results. Ten liters of 1x buffer can be made from 10x buffer stock.

Premixed buffers are made from the high quality reagents and are carefully mixed and tested to ensure reproducibility.

10x Composition

0.25 M Tris base
1.92 M Glycine
pH approx. 8.3

Recommended Storage Condition

Store at room temperature

Usage Information (applications)

Electrophoresis buffer for Tris-Glycine gel and western blotting

Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.



For research use only

ISO9001

10x Tris-Glycine Running Buffer

Product	Volume	Cat. No.
10x Tris-Glycine Running Buffer	1 L (10x)	EBA-1042

Description.

Elpis Biotech's Tris-Glycine Running Buffer (10x) is suitable for use as an electrophoresis buffer for Tris-Glycine gel to separate proteins and for western blotting. Premixed electrophoresis buffer eliminate variables which can distort data and give unreliable results. Ten liters of 1x buffer can be made from 10x buffer stock.

Premixed buffers are made from the high quality reagents and are carefully mixed and tested to ensure reproducibility.

10x Composition

0.25 M Tris base
1.92 M Glycine
pH approx. 8.3

Recommended Storage Condition

Store at room temperature

Usage Information (applications)

Electrophoresis buffer for Tris-Glycine gel and western blotting

Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.