



For research use only

ISO9001

1M Tris-HCl (pH7.0, 7.5, 8.0)

Product	Volume	Cat. No.
1 M Tris-HCl (pH7.0)	500 ml	EBA-1102
(pH7.5)	500 ml	EBA-1103
(pH8.0)	500 ml	EBA-1104

Description.

Elpis Biotech's 1 M Tris-HCl (pH7.0, 7.5, 8.0) is made from the high quality Tris base and are filtered by 0.2 μ m filter. Tris-HCl, Molecular Biology Grade, is sometimes used in combination with Tris base for preparation of Tris-HCl buffers

Composition

1 M Tris base
pH 7.0, 7.5, 8.0 at 25°C

Recommended Storage Condition

Store at 2-8°C

Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.



For research use only

ISO9001

1M Tris-HCl (pH7.0, 7.5, 8.0)

Product	Volume	Cat. No.
1 M Tris-HCl (pH7.0)	500 ml	EBA-1102
(pH7.5)	500 ml	EBA-1103
(pH8.0)	500 ml	EBA-1104

Description.

Elpis Biotech's 1 M Tris-HCl (pH7.0, 7.5, 8.0) is made from the high quality Tris base and are filtered by 0.2 μ m filter. Tris-HCl, Molecular Biology Grade, is sometimes used in combination with Tris base for preparation of Tris-HCl buffers

Composition

1 M Tris base
pH 7.0, 7.5, 8.0 at 25°C

Recommended Storage Condition

Store at 2-8°C

Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.