



For research use only

ISO9001

## Nucleic Acid Purification Column

Product	Qty	Cat. No.	Remarks
Spin Column (mini)	200 ea	EBD-1021	silica-filter in spin column

### Description

Elpis Biotech's Nucleic Acid Purification columns are suitable for purification of plasmid DNA, PCR product and genomic DNA.

### Binding capacity

Spin Column	200 ea with collection tubes (2 ml)
Binding Capacity	40-50 $\mu$ g

### Recommended Storage Condition

Store at Room Temperature.

### Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.



For research use only

ISO9001

## Nucleic Acid Purification Column

Product	Qty	Cat. No.	Remarks
Spin Column (mini)	200 ea	EBD-1021	silica-filter in spin column

### Description

Elpis Biotech's Nucleic Acid Purification columns are suitable for purification of plasmid DNA, PCR product and genomic DNA.

### Binding capacity

Spin Column	200 ea with collection tubes (2 ml)
Binding Capacity	40-50 $\mu$ g

### Recommended Storage Condition

Store at Room Temperature.

### Caution

Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.