



For research use only

ISO9001

Afl II

Product	Quantity	Cat. No.	Remarks
Afl II	2,500 unit	EBR-1201	10 unit/ μ l

Recognition site**Source***Anabaena flos-aquae*.**Reaction Conditions**Buffer 4 + BSA. 20 mM Tris-acetate (pH 7.9), 50 mM potassium acetate, 10 mM magnesium acetate, 1 mM DTT. Supplement with 100 μ g/ml BSA. Incubate at 37°C.**Concentration & Storage Condition**Store at -20°C. 10 unit/ μ l in 10 mM Tris-HCl (pH 7.4), 50 mM KCl, 1 mM DTT, 0.5 mM EDTA, 200 μ g/ml BSA, 50% (v/v) glycerol.**Heat Inactivation Condition**

65°C for 20min.

Enzyme dilution

Dilute with Storage Buffer.

QC Tests

Activity, ligation and recutting, exo and endonuclease activity test, SDS-PAGE purity, performance tests.

Methylase SensitivityNot sensitive to *dam*, *dcm* or CpG methylation.

Buf4	37°C	65°C		BW
<small>BUFFER</small>	<small>REACTION</small>	<small>INACTIVATION</small>	<small>METHYL BLOCK</small>	<small>BLUE/WHITE</small>

Buffer	1	2	3	4
Activity	50	100	25	100



For research use only

ISO9001

Afl II

Product	Quantity	Cat. No.	Remarks
Afl II	2,500 unit	EBR-1201	10 unit/ μ l

Recognition site**Source***Anabaena flos-aquae*.**Reaction Conditions**Buffer 4 + BSA. 20 mM Tris-acetate (pH 7.9), 50 mM potassium acetate, 10 mM magnesium acetate, 1 mM DTT. Supplement with 100 μ g/ml BSA. Incubate at 37°C.**Concentration & Storage Condition**Store at -20°C. 10 unit/ μ l in 10 mM Tris-HCl (pH 7.4), 50 mM KCl, 1 mM DTT, 0.5 mM EDTA, 200 μ g/ml BSA, 50% (v/v) glycerol.**Heat Inactivation Condition**

65°C for 20min.

Enzyme dilution

Dilute with Storage Buffer.

QC Tests

Activity, ligation and recutting, exo and endonuclease activity test, SDS-PAGE purity, performance tests.

Methylase SensitivityNot sensitive to *dam*, *dcm* or CpG methylation.

Buf4	37°C	65°C		BW
<small>BUFFER</small>	<small>REACTION</small>	<small>INACTIVATION</small>	<small>METHYL BLOCK</small>	<small>BLUE/WHITE</small>

Buffer	1	2	3	4
Activity	50	100	25	100